



**BRIZ 2** (Code:SW-8021P)

- Steel powder coating
- Seat width : 16" / 18"
- Weight capacity : 100 kg
- Detachable armrest
- Detachable footrest
- 8" solid caster
- 24" solid rear wheel

Optional :

- Foldable backrest (Briz 3)

**BRIZ 3** (Code:SW-8520P)



Foldable backrest



Model									
SW-8021P	41/46cm	56/64cm	27cm	106W	50cm	40cm	95cm	17	100
SW-8520P	41/46cm	56/64cm	27cm	106W	50cm	40cm	95cm	17	100

±1cm ± 0.5kg

- KARMA reserves the right to modify information herein without further notice.