

HAND 701

KM-7501



Ultra lightweight and compact



Here-we-go handle



Side-lifting transparent tray



Features

- ✓ Ergonomically antelope horn-shaped handle makes it easy to steer and push the chair
- ✓ Ergonomically designed vertical footrest allow legs to be placed in the correct position
- ✓ Outward extended front wheels
- ✓ One-piece footplate increases stability ultra lightweight and compact

Specifications

- ✓ Seat width : 11" / 13.5"
- ✓ Weight capacity : 60 kg
- ✓ Fixed armrest ✓ Fixed footrest
- ✓ 6" solid caster
- ✓ 14" solid rear wheel

Width	F/R Wheels	Seat Depth	Seat Width	Seat Height	Backrest Height	Backrest Angle	Wheelchair Weight	Wheelchair Height	Wheelchair Width	Wheelchair Length	Wheelchair Depth	Max. User Weight
11"	6"-14F"	28	30	45	34	18	70	39	36	102	9.3	60kg
13.5"	6"-14F"	28	30	45	34	18	70	39	36	102	9.3	60kg

• KARMA reserves the right to modify information herein without further notice. ±1cm ± 0.5kg