

Aluminum Alloy Wheelchair

KM-8020X



Features

- ✓ 160 Kg weight loading capacity.
- ✓ Available in 20" & 22" seat width.
- ✓ Crash tested to ISO 7179-1917.
- ✓ Foldable backrest & swing-away footrest.

Specifications

- ✓ Seat width : 20" / 22"
- ✓ Weight capacity : 160 kg
- ✓ Foldable backrest
- ✓ Flipback armrest
- ✓ Wide profile casters
- ✓ CE certified

Width	F/R Wheels										
22"	8"-24"	56.5cm	46.5cm	76.5cm	37cm	23cm	105cm	53cm	40cm	93cm	16.9Kg, 160kg (350lbs)
<small>• KARMA reserves the right to modify information herein without further notice. ±1cm ± 0.5kg</small>											



KARMA[®]

Mobility Solution

Dear Sir,

We at Karma Healthcare India have always strived to deliver our range of products which are designed as per the World Class Standards for which KARMA as a brand is known for.

We have immense pleasure in introducing to you our new arrival KM-8020X, its now available for order booking in India.

The Key features of KM-8020X are :

- 160 Kg weight loading capacity.
- Available in 20' & 22' seat width.
- Crash tested to ISO 7179-1917.
- Foldable backrest & swing-away footrest.

Our Sales Team would get in touch with you for trade prices and other details.

Once again, as we always say at KARMA – Discover the path and make it yours!

Wish you all the very best in business and life!!
Assuring you of our best attention at all times.

Sincerely yours
For KARMA HEALTHCARE INDIA

D.N.Trivedi
Chief Executive Officer
Karma Healthcare India Ltd

KARMA HEALTHCARE LIMITED

INDIA

GD 347, Block-GD, Sector-3, Salt Lake City, Kolkata - 700 106, West Bengal, India

Phone - 033 4022 1900 ; Fax - 033 4022 1910 ; email - sganguly@karmahealthcare.com ; web - www.karmahealthcare.com