

# Power indoor Wheelchair

(Code:KP-10.2)



Option: Curb climber



Quick attachable/  
detachable battery  
design



Super compact  
to fit your trunk



Footrests can be  
moved forward or  
backward for optimum  
placement



Height adjustable  
armrest



Adjustable armrest  
width



KP-10.2 can be disassembled  
easily for transportation and  
storage.

## Colors



Silvery Black

Width																
KP-10.2	90cm	62cm	97cm	43cm	F:7.8x2" R:12.5x2.25"	180W	22Ahx2	Shark/VR2	76cm	12"	6km	16-20km	Aluminum Alloy	115kg	49.46kg	34.76kg
± 1cm ± 0.5kg																

The above are estimated data under 20°C (68°F), an 80kg (176lbs) user, brand new fully charged batteries and a constant driving speed of 6km/hr (3.7mph). Distance, climbing ability and operating time can vary according to condition.

KARMA reserve the right to modify information herein without further notice.

## KP-10.2 Features

- ⊙ Seat, base & battery quickly detach for easy storage and transportation
- ⊙ With seat platform, the captain seat can be moved forward or backward using tools included.
- ⊙ Adjustable armrest width and height
- ⊙ Footrests can be moved forward or backward for optimum placement
- ⊙ Quick release battery pack for easy storage and battery pack change
- ⊙ Controller can be easily switched for left or right hand usage
- ⊙ The heaviest component is merely 24 kg (base)
- ⊙ Interchanges with standard KP-10 seat

## Options

- ⊙ Upper cushion for back & seat upholstery (KP-10)
- ⊙ Curb Climber (KP-10.2)