

# KP 25.2



59.5cm	104cm	40.6cm	F:7.8x2" R:12.5x2.25"	180Wx2	36Ahx2	VR2	N/A	95.5cm	90cm	12"	6km	30Km	Aluminum Alloy	115kg	58.6kg
± 1cm ± 0.5kg															

The above are estimated data under 20°C (68°F), an 80kg (176lbs) user, brand new fully charged batteries and a constant driving speed of 6km/hr (3.7mph). Distance, climbing ability and operating time can vary according to condition.

KARMA reserve the right to modify information herein without further notice.

- Features**
- ⊙ Aluminum-alloy frame that won't rust
  - ⊙ Easy access battery carrying system
  - ⊙ Foldable backrest to allow compact transportation and storage
  - ⊙ Built-in multi-functional controller
  - ⊙ Ultra-lightweight design makes transport easy and quick
  - ⊙ The heaviest component is merely 32 kg

- Options**
- ⊙ Curb climber
  - ⊙ Foldable drink holder
  - ⊙ Attendant controller system

