

Power Wheelchair

KP 45.3 TR



Features

- ✓ Power Tilt & Recline relieve pressure
- ✓ Model of choice for active outdoor use
- ✓ Multi function controller
- ✓ Quick release seat for easy maintenance
- ✓ Comfortable sitting position reduce the risk of decubitus ulcer

Specifications

- ✓ 400 W high performance motor
- ✓ Weight capacity : 135 kg
- ✓ Total weight 107 kg with battery
- ✓ Ergonomically designed seats
- ✓ Suspension and anti tippers



18"20



65cm



127cm



122cm



52cm



400wx
2 pcs



76 kg



107 kg



135

KARMA reserves the right to modify information herein without further notice.

±1cm ± 0.5kg

Unit : cm



Dear Sir,

We at Karma Healthcare India have always strived to deliver our range of products which are designed as per the World Class Standards for which KARMA as a brand is known for.

We have immense pleasure in introducing to you our new arrival KP 45.3 TR, its now available for order booking in India.

The Key features of KP 45.3 TR are :

- Ergonomically designed seat.
- 400 W high performance Motor.
- Suspension anti-tippers.
- Multi-functional controller.
- Power Tilt & Recline relieve pressure (Multi-functional seat)
- Comfortable sitting position reduce the risk of decubitus ulcer

Our Sales Team would get in touch with you for trade prices and other details.

Once again, as we always say at KARMA – Discover the path and make it yours!

Wish you all the very best in business and life!!
Assuring you of our best attention at all times.

Sincerely yours
For KARMA HEALTHCARE INDIA

D.N.Trivedi
Chief Executive Officer
Karma Healthcare India Ltd

KARMA HEALTHCARE LIMITED

INDIA
GD 347, Block-GD, Sector-3, Salt Lake City, Kolkata - 700 106, West Bengal, India
Phone - 033 4022 1900 ; Fax - 033 4022 1910 ; email - sganguly@karmahealthcare.com ; web - www.karmahealthcare.com