


■ Fully adjustable seat & armrest system for great comfort on your demand



Height Adjustable



Forward/Backward Adjustable



Width Adjustable



Detachable



Freewheel hub



Suspension anti-tippers



Built-in lighting system



Quick-release seat design



Multi-functional controller



Hand bar(Optional)



Push bar(Optional)

118cm	65cm	118cm	50cm (20")	400Wx2	50Ahx2	Shark60A	Rear Wheel/ Anti-tipper	F: 8.6" R: 16"	90cm	12'	13km/hr	40-45km	135kg	99.5kg	68kg
														± 1cm ± 0.5kg	

The above are estimated data under 20°C (68°F), an 80kg (176lbs) user, brand new fully charged batteries and a constant driving speed of 10km/hr (6.2mph). Distance, climbing ability and operating time can vary according to condition.

KARMA reserve the right to modify information herein without further notice.

• Features

- Quick-release seat design provides easy maintenance
- Unique 2-D suspension system
- Suspension anti-tippers
- Built-in lighting system
- Adjustable seat & armrest system
- Additional built in charger socket to protect controller
- Recline backrest

Options

- High-speed tires
- Elevating footrest
- Parallel swing-away Joystick bracket
- 75Ah batteries
- Curb climber (please notify us upon placing order)
- Drink holder
- Hand bar
- Push bar

Colors

