

## KS -343



(4-wheel model)



Option: Double-pivot steering column



Taillight, turn signals & bumper

Width															
KS-343	62cm	112cm	116cm	350W	36Ahx2	5Amp	F/R	F/R: 10"	120cm	12□	9.3km/hr	35km	135kg	92.2kg	66.8kg
KS-333	62cm	112cm	116cm	350W	36Ahx2	5Amp	F/R	F/R: 10"	110cm	12□	9.3km/hr	35km	135kg	90kg	66kg

The above are estimated data under 207(68?), an 80kg(176lbs) user, brand new fully charged batteries and a constant driving speed of 9.3km/hr (5.7mph). Distance, climbing ability and operating time can vary according to condition. ± 1cm ±0.5kg

## The best partner for around town or in the mall

The KS-343 is an economically and ingeniously designed compact scooter, whose chic abounds.

Its ergonomically designed footrest gives more than enough legroom for taller individuals.

The rear hydraulic suspension system provides impressive stability and a smooth ride. It's easy to disassemble the front and rear frames with its quick-release lever.

The KS-343 is an excellent choice for driving in the city, around the neighborhood or in a shopping mall.

### Features

- ⊙ Energy saving LED taillights
- ⊙ Stylish aluminum-alloy wheels
- ⊙ Larger rear wheels for added stability
- ⊙ Rear shock absorber system
- ⊙ Front bumper for added safety
- ⊙ Double-pivot steering column for option

### Colors



Ruby red



Silver pearl



Ocean green