

# KS-747.2



■ Reflective film



■ Rear shock absorber & bumper



■ Taillight and turn signals



■ Double-pivot steering column (option)



■ Headlight and front bumper



■ Front shock absorber

Width															
KS-747.2	64.5cm	127.5cm	110-117cm	550W	50Ahx2	8Amp	F/R	F/R:12"	166cm	15"	15km/hr	50km	135kg	105.9kg	74.3kg

The above are estimated data under 20°C (68°F), an 80kg (176lbs) user, brand new fully charged batteries and a constant driving speed of 10km/hr (6.2mph). Distance, climbing ability and operating time can vary according to condition.

± 1cm ± 0.5kg

KARMA reserve the right to modify information herein without further notice.

## Colors



Ruby red



Silver pearl



Ocean green

## Features

- Front and rear shock absorber systems provide user great comfort
- Double-pivot steering column optional
- The wide-angle headlight on the nameplate casing helps lead the way by maintaining its angle when adjusting the steering column
- Well-ventilated synthetic leather seat with adjustable height for greater comfort
- Taillight and front & rear turn signals enhance driving safety