

# Standard Steel Wheelchair



## **SUNNY 1** Powder Coating frame (Code:SW-1000P)

- Seat width : 16" / 18"
- Weight capacity : 120 kg
- Fixed armrest
- Fixed footrest
- 8" solid caster
- 24" solid rear wheel

## **SUNNY 1**

Chromed Plated frame (Code:SW-1000C)



Model									
SW-1000P	46/41cm	60/65cm	27cm	104W	47cm	42cm	91cm	21	120
SW-1000C	46/41cm	60/65cm	27cm	104W	47cm	42cm	91cm	21	120

\* KARMA reserves the right to modify information herein without further notice. ±1cm ± 0.5kg

## **SUNNY 3** Powder Coating frame (Code:SW-8040P)

- Seat width : 16" / 18"
- Weight capacity : 100 kg
- **Detachable armrest**
- **Detachable footrest**
- 8" solid caster
- 24" solid rear wheel
- Elevating footrest



Elevating footrest